



## MICRO-HEM82

### Fully Auto 3-Part Hematology Analyzer

Hematology Analyzer Fully automatic

9.7-inch color touch LCD screen, easy to operate, display all parameters and histograms on the same screen.

Fewer samples are required for testing: whole bloods 12µl, pre-diluteds 20µl.

60 samples/hour.

21 parameters + 3 histograms.

It can be turned on continuously for 24 hours.

Automatically clean the sampler, counting pool and pipeline after starting up.

Electrical impedance method, cyanide-free detection of HGB

Both automatic calibration and manual calibration are available.

Counting Mode: Venous, capillary, prediluted. @. Support LIS system.

- **On request veterinary software available**



## MICRO-HEM85

### Fully Auto 5-Part Hematology Analyzer

5 diff differential of WBC

Double chambers for WBC and RBC/ PLT

27 parameters

4 histograms

1 scatter diagram

Up to 90 samples/hour

Touch screen display

Storage: up to 100 thousand results with graphs

Barcode scanner for option

Dimensions & Weight: 53.2x 40.0x 52.0 cm, 28 kgs

- **On request Veterinary software available**
- **On request with CRP**
- **On request with Autoloader**